

**IR GM300 Thermomtre infrarouge missivit fixe (0,95)
(usage non mdical uniquement)**



[Click here to view product](#)

Product Code	XE-4008-001
Model	GM300 Fixed Emissivity (0.95)
General Description	GM300 infrared thermometer is used for measuring the temperature of the object's surface.
Accuracy	+/- 1.5Å°C / +/- 1.5%
Resolution	0.1Å°C or 0.1Å°F
Emissivity	0.95 Preset
Distance to Spot Size	12:1
Battery	9V
Function #1	MAX / MIN / AVG / DIF reading



Function #2	Auto Power off & Data Hold
Function #3	High / Low Temperature Alarm setup
Function #4	Laser ON / OFF selectable
Unit selectable	°C / °F
Power Supply	Requires one standard 9V battery (included)
Measuring range	-50°C to 450°C
Important Information	These are industrial IR thermometers, so please check the specifications before purchasing as a restocking fee of 20% will apply to unused products