


## Sonde Pt100 usage gnral avec fil PFA - Type RTD / PRT



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Sonde Pt100 usage gnral avec fil PFA - Type RTD / PRT**

Sonde polyvalente Pt100 selon IEC751

Classe B, connexion 4 fils

Gaine rigide en acier inoxydable 316, 3,0, 4,0, 4,5 ou 6,0 mm de diamtre

Plomb blind isol au tflon de 1 ou 2 mtres

Temprature nominale du capteur -50C +250C

**Spcifications****Specifications**

<b>General Description</b>	Multipurpose Pt100 Probe to IEC751
<b>Sensor Type</b>	Pt100 (100 Ohms @ 0°C) to IEC 751, Class B, 4 wire
<b>Sheath Material</b>	3.0, 4.0, 4.5 or 6.0mm diameter 316 stainless steel sheath.
<b>Cable Length</b>	1 or 2 metre flexible Teflon® insulated, screened
<b>Termination Type</b>	4-wire tails allows connection to any PT100 instrument 2, 3 or 4 wire
<b>Max. Temperature</b>	+250°C
<b>Min. Temperature</b>	-50°C
<b>Standards Met</b>	IEC751