


Sonde Pt100 usage gnral avec fil PFA - Type RTD / PRT



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Sonde Pt100 usage gnral avec fil PFA - Type RTD / PRT

Sonde polyvalente Pt100 selon IEC751

Classe B, connexion 4 fils

Gaine rigide en acier inoxydable 316, 3,0, 4,0, 4,5 ou 6,0 mm de diamtre

Plomb blind isol au tflon de 1 ou 2 mtres

Température nominale du capteur -50C +250C

Spcifications**Specifications**

General Description	Multipurpose Pt100 Probe to IEC751
Sensor Type	Pt100 (100 Ohms @ 0°C) to IEC 751, Class B, 4 wire
Sheath Material	3.0, 4.0, 4.5 or 6.0mm diameter 316 stainless steel sheath.
Cable Length	1 or 2 metre flexible Teflon® insulated, screened
Termination Type	4-wire tails allows connection to any PT100 instrument 2, 3 or 4 wire
Max. Temperature	+250°C
Min. Temperature	-50°C
Standards Met	IEC751