



# LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

## Raccords compression en acier inoxydable - Filetage mtrique (srie M)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Pour adapter une large gamme de capteurs isolation minérale et autres capteurs de type sonde proposés par Labfacility.

Les raccords compression en acier inoxydable filetage métrique (série M) sont largement utilisés pour le raccordement de tuyaux et de tubes dans diverses applications industrielles.

Ces raccords sont dotés de filetages métriques, normalisés pour assurer la compatibilité avec les systèmes de tuyauterie métriques.

Fabriqués en acier inoxydable, ils offrent une résistance à la corrosion, une durabilité et une résistance supérieures, ce qui les rend adaptés aux environnements difficiles tels que le traitement chimique, les produits pharmaceutiques et les industries agroalimentaires.

### Spécifications

## Specifications

<b>General Description</b>	Thermocouple installation fittings with Metric Thread
<b>Material</b>	Male in-line stainless steel compression type with nut & olive
<b>Thread</b>	M8, M16 & M20 thread available.
<b>Probe Diameter</b>	Suits 1.0mm, 1.5mm, 2.0mm, 3.0mm, 4.5mm, 6.0mm or 8.0mm probe diameters