

## Sonde de serrage de tube - type K



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Cette sonde de serrage convient la mesure de la temprature de surface des tuyaux dans les systmes de chauffage et de ventilation, la rfrigration et de nombreuses autres applications de tuyauterie.

## Specifications

<b>Product Code</b>	XE-3847-001
<b>General Description</b>	Pipe Clamp Probe
<b>Sensor Type</b>	Clamps onto pipes from 6 mm to 30 mm
<b>Thermocouple Type</b>	K
<b>Cable Length</b>	1m
<b>Cable Type</b>	PVC
<b>Cable Termination</b>	Miniature Plug
<b>Max. Temperature</b>	-10°C
<b>Min. Temperature</b>	+100°C
<b>Accuracy</b>	Response time less than 2 seconds