


## Platinum Resistance Pt100 Sonde industrielle avec extension de lagging



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Besoin de votre capteur de température calibré?**

Traçable aux normes d'étalonnage ISO17025 (UKAS), Labfacility offre un choix d'un étalonnage traçable de 3 ou 5 points à l'aide de notre installation d'étalonnage interne, vous pouvez sélectionner tous les points de température entre -10 °C et 200 °C.

**Veillez noter que si vous ajoutez une option d'étalonnage à votre commande, vous serez contacté concernant les informations de livraison.**

**Specifications**

<b>Sensor Type</b>	Pt100 (100 Ohms @ 0°C) to IEC 751
<b>Insert</b>	6.0mm diameter 316 stainless steel terminated with spring loaded ceramic connection block
<b>Process Connection</b>	1/2"BSP parallel
<b>Cable Connection</b>	M20 x 1.5mm cable entry (gland included)
<b>Termination Type</b>	IP67 aluminium alloy weatherproof KNE connection head
<b>Max. Temperature</b>	+450°C (connection head @ 170°C)
<b>Min. Temperature</b>	-50°C
<b>Tolerance</b>	Class B
<b>Number of Wires</b>	4
<b>Probe Material</b>	8.0mmO 316 stainless steel below 1/2"BSP parallel process connection 1/2" diameter x 75mm lagging length to underside of connection head