


Manomètre de pression GM520



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

[Cliquez ici](#)

pour télécharger les derniers logiciels nécessaires a l'exploitation de ce produit

Spécifications**Specifications**

Product Code	XE-4011-001
Model	GM520 Pressure Manometer
Communications Protocol	USB communication
Accuracy	+/- 0.3% FSO (25°C)
Repeatability	+/- 0.2% (Maximum +/- 0.5% FSO)
Response	Typical 0.5 s
Linearity/hysteresis	+/- 0.29% FSO
Maximum Pressure	+/- 150kPa
Storage Temperature Range	-10 to 60°C
Ambient Operating Temperature	0 to 50°C
Measuring range	+/- 35kPa