


## Manomètre de pression GM520



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

[Cliquez ici](#)

pour télécharger les derniers logiciels nécessaires a l'exploitation de ce produit

**Spécifications****Specifications**

<b>Product Code</b>	XE-4011-001
<b>Model</b>	GM520 Pressure Manometer
<b>Communications Protocol</b>	USB communication
<b>Accuracy</b>	+/- 0.3% FSO (25°C)
<b>Repeatability</b>	+/- 0.2% (Maximum +/- 0.5% FSO)
<b>Response</b>	Typical 0.5 s
<b>Linearity/hysteresis</b>	+/- 0.29% FSO
<b>Maximum Pressure</b>	+/- 150kPa
<b>Storage Temperature Range</b>	-10 to 60°C
<b>Ambient Operating Temperature</b>	0 to 50°C
<b>Measuring range</b>	+/- 35kPa