

GM1361 Compteur de temprature et dhumidit



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Actuellement en attente de nouveau stock - pr-commande uniquement.

Vous serez averti lorsque votre commande sera prête à être expédiée.

Specifications

Product Code	XE-4012-001
Model	GM1361 Temperature & Humidity Meter
General Description	Supplied in hard carry case
Type	K-type thermocouple temperature measuring
Accuracy of temperature	±1°C or ±1.8°C
Accuracy of humidity	± 3%(30-95%); ± 5%(10-30%)
Battery	9V Battery
Dimensions	160*70*30mm
Function #1	Large LCD temperature & humidity dual display
Function #2	Data hold
Function #3	Low battery indication
Function #4	Auto /Manual power shut off
Function #5	Max/ Min function
Input Range - Type K	-10°C ~ 50°C (-30°C ~ 1000 °C K-type)
Humidity range	5%RH ~ 98%RH
Unit selectable	°C / °F selection
Sampling rate	2.5times/sec