





# LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

## Sonde de température de pneus de course réglable



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

## Specifications

<b>Product Code</b>	XE-3892-001
<b>General Description</b>	This very fast response racing tyre probe has a manual depth adjustment from 1 mm to 15 mm
<b>Sensor Type</b>	Can be used with any type K digital thermometer with a mini input socket.
<b>Thermocouple Type</b>	K
<b>Probe Diameter</b>	1.1mm
<b>Probe Length</b>	manual depth adjustment from 1 mm to 15 mm
<b>Cable Length</b>	1m
<b>Cable Termination</b>	Miniature Plug
<b>Max. Temperature</b>	+250°C
<b>Min. Temperature</b>	-50°C
<b>Features #1</b>	Junction Type: Grounded