




LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Sonde de température de pneus de course réglable



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications

Product Code	XE-3892-001
General Description	This very fast response racing tyre probe has a manual depth adjustment from 1 mm to 15 mm
Sensor Type	Can be used with any type K digital thermometer with a mini input socket.
Thermocouple Type	K
Probe Diameter	1.1mm
Probe Length	manual depth adjustment from 1 mm to 15 mm
Cable Length	1m
Cable Termination	Miniature Plug
Max. Temperature	+250°C
Min. Temperature	-50°C
Features #1	Junction Type: Grounded