

**4-20mA boÃfÃ®tier ÃfÃ distance montÃfÃ© par
mur, entrÃfÃ©e Pt100 avec PFA ou cÃfÃble
SSOB - capteur PRT**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Logement moulé dans l'ABS a usage général, les points de fixation externes fournissent un moyen simple et rapide de montage des capteurs directement sur les conduits, dispensant de la nécessité de plaques de montage séparées.

Besoin de votre capteur de température calibré?

Traçable aux normes d'étalonnage ISO17025 (UKAS), Labfacility offre un choix d'un étalonnage traçable de 3 ou 5 points a l'aide de notre installation d'étalonnage interne, vous pouvez sélectionner tous les points de température entre -10 oC et 200 oC.

Veillez noter que si vous ajoutez une option d'étalonnage a votre commande, vous serez contacté concernant les informations de livraison.

Specifications

General Description	Compatible with Pt100 platinum resistance thermometers
Sensor Type	Surface Mounting for remote installation
Cable Length	1 or 2 metre
Cable Type	PFA screened or Fiberglass stainless steel overbraid available
Core / Strands	7/0.2mm
Head Type	High Impact ABS Plastic Construction
Thread	Entry Threads M16
Transmitter	Yes
Range	-50°C to +150°C or 0°C to +400°C (either can be re-ranged if desired)
Output	4-20 mA output