

**Botier mural distance 4-20mA, entre Pt100 avec
cble PFA et capteur PRT - Longueur de cble 1M**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Botier moule en ABS usage gnral, les points de fixation externes offrent un moyen simple et rapide de monter des capteurs directement sur les conduits, sans avoir besoin de plaques de montage spares.

Vous avez besoin de calibrer votre capteur de temprature ?

Traable selon les normes dtalonnage ISO17025 (UKAS), Labfacility offre un choix dtalonnage traable en 3 ou 5 points laide de notre installation dtalonnage interne, vous pouvez slectionner nimporte quel point de temprature entre -10C et +200C.

Veillez noter que si vous ajoutez une option de calibrage votre commande, vous serez contact pour les informations de livraison.

Spifications

Specifications

Product Code	XE-5660-001
General Description	Compatible with Pt100 platinum resistance thermometers
Sensor Type	Surface Mounting for remote installation
Probe Diameter	6mm
Probe Length	250mm
Cable Length	1m
Cable Type	PFA screened
Core / Strands	7/0.2mm
Head Type	High Impact ABS Plastic Construction
Thread	Entry Threads M16
Transmitter	Yes
Range	-50°C to +150°C
Output	4-20 mA output