

## Botier mural distance 4-20mA, entre Pt100 avec cble SSOB et capteur PRT - longueur de cble 2M



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Botier moul en ABS usage gnral, les points de fixation externes offrent un moyen simple et rapide de monter des capteurs directement sur les conduits, sans avoir besoin de plaques de montage spares.

### **Vous avez besoin de calibrer votre capteur de temprature ?**

Traable selon les normes dtalonnage ISO17025 (UKAS), Labfacility offre un choix dtalonnage traable en 3 ou 5 points laide de notre installation dtalonnage interne, vous pouvez slectionner nimporte quel point de temprature entre -10C et +200C.

**Veillez noter que si vous ajoutez une option de calibrage votre commande, vous serez contact pour les informations de livraison.**

### **Spcifications**

## **Specifications**

|                            |  |
|----------------------------|--|
| <b>Product Code</b>        | XE-5663-001  |
| <b>General Description</b> | Compatible with Pt100 platinum resistance thermometers |
| <b>Sensor Type</b>         | Surface Mounting for remote installation               |
| <b>Probe Diameter</b>      | 6mm  |
| <b>Probe Length</b>        | 300mm  |
| <b>Cable Length</b>        | 2m   |
| <b>Cable Type</b>          | Fiberglass stainless steel overbraid available         |
| <b>Core / Strands</b>      | 7/0.2mm  |
| <b>Head Type</b>           | High Impact ABS Plastic Construction                   |
| <b>Thread</b>              | Entry Threads M16                                      |
| <b>Transmitter</b>         | Yes  |
| <b>Range</b>               | 0°C to +400°C  |
| <b>Output</b>              | 4-20 mA output   |