

## Thermomtre numrique de restauration type T



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Thermomtre de restauration**

Thermomtre de chef de haute prcision

Thermocouple type T

Conforme HACCP,

-100 + 280C

Fonction de maintien,

Thermomtre tanche IP67

Double cran LCD

Facile utiliser: fonctionnement deux boutons

Double affichage

Capacit darrt automatique.

Fonctionnalits:

Auto-talonnage

Indication de batterie faible

Livr avec 1 batterie PP3 9V

Plage de mesure -100 +280C

Haute prcision +/- 0,1 % de la lecture

+/-0,2C, 0,1C Rsolution

Le CA2005 offre une fonction de maintien et lorsque vous appuyez sur ce bouton, la temprature actuelle sera maintenue et transfre en haut droite de lcran. Lcran principal continuera dafficher la temprature en temps rel.

**Spfications****Specifications**

<b>Product Code</b>	XE-4090-001
<b>General Description</b>	High accuracy chef thermometer
<b>Type</b>	IP67 waterproof thermometer



# LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

<b>Accuracy</b>	+/-0.1% of Reading
<b>Resolution</b>	+/-0.2°C, 0.1°C Resolution
<b>Input / Ranges</b>	-100 to + 280°C
<b>Power Supply</b>	Supplied with 1 9V PP3 battery
<b>Sensor Type</b>	T