

Cble thermocouple monoshot PTFE Type K



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

PTFE Single Shot Thermocouple Type K

Construction unique de cble monoshot
Rsisant aux produits chimiques, aux huiles et aux fluides. Convient aux applications arospatiales et est largement
utilis dans les autoclaves.

Spifications

Specifications

General Description	Unique single shot cable construction Resistant to chemicals, oils and fluids. Suitable for aerospace applications and is widely used in autoclaves.
Cable Length	10, 25, 50, 100M
Cable Type	PTFE
Core / Strands	1/0.376mm
Colour	Green
Max. Temperature	+250°C
Standards Met	IEC-584-3
Thermocouple Type	K