

Cble thermocouple monoshot PTFE Type T



50M



10m



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

PTFE Single Shot Thermocouple Type T

Construction unique de cble monoshot
Rsisant aux produits chimiques, aux huiles et aux fluides. Convient aux applications arospatiales et est largement
utilis dans les autoclaves.

Specifications

Specifications

| | |
|----------------------------|--|
| General Description | Unique single shot cable construction Resistant to chemicals, oils and fluids. Suitable for aerospace applications and is widely used in autoclaves. |
| Cable Length | 10, 25, 50, 100M |
| Cable Type | PTFE |
| Core / Strands | 1/0.376mm |
| Colour | Brown |
| Max. Temperature | +250°C |
| Standards Met | IEC-584-3, Class 1 |
| Thermocouple Type | T |