

## Tte de raccord KNE



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

De nombreux types alternatifs de ttes de raccordement sont disponibles pour rpondre aux exigences de diverses applications.

Il existe des variations en termes de taille, de matriau, dhbergement, de rsistance aux fluides, de rsistance au feu ou mme lexplosion et dautres paramtres.

### Spfications

## Specifications

<b>Head Type</b>	KNE Standard Head
<b>Material</b>	Aluminium silver painted
<b>Cable Entry</b>	M20 x 1.5mm
<b>Fixing (Centres)</b>	33mm & 40mm
<b>Max. Temperature</b>	170°C
<b>IP Rating</b>	IP67 rated
<b>Process Entry</b>	1/2" BSPP, 3/4" BSPP or M24
<b>Cap Retaining Type</b>	Screw cap and retaining screw