

Tte de raccord KNE



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

De nombreux types alternatifs de ttes de raccordement sont disponibles pour rpondre aux exigences de diverses applications.

Il existe des variations en termes de taille, de matriau, dhbergement, de rsistance aux fluides, de rsistance au feu ou mme lexplosion et dautres paramtres.

Spfications

Specifications

Head Type	KNE Standard Head
Material	Aluminium silver painted
Cable Entry	M20 x 1.5mm
Fixing (Centres)	33mm & 40mm
Max. Temperature	170°C
IP Rating	IP67 rated
Process Entry	1/2" BSPP, 3/4" BSPP or M24
Cap Retaining Type	Screw cap and retaining screw