

**Capteurs Pt1000 Classe B avec cble isol en tflon -
Type PRT - 4,0 mm de diamtre - 40 mm de
longueur**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications**Specifications**

Product Code	XE-3682-001
General Description	Insulated element in a 4.0mm diameter 316 stainless steel probe
Sensor Type	Pt1000 (1000 Ohms @ 0°C) to IEC 751, Class B
Element	Pt1000
Probe Diameter	4.0mm
Probe Length	40mm
Sheath Material	316 stainless steel probe
Cable Length	1M
Cable Type	Teflon® insulated twisted leads
Cable Termination	Bare tails
Max. Temperature	200°C
Min. Temperature	-50°C
Number of Wires	2
Series	General Purpose Pt1000