


Type J Glassfibre Insulated Flat Pair thermocouple Cable / Fil (BS)

Type J, 1/0.508mm
Fiberglass Flat pair



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Type J Glassfibre Insulated Flat Pair thermocouple Câble / Fil (BS)

L'image est pour l'illustration de codage de couleur seulement

Le prix est par metre - s'il vous plaît ajouter la quantité nécessaire pour le panier. Il ne reste plus que 1 227 metres en stock

Construction isolée de paires plates de Glassfibre, chefs d'orchestre posés a plat, Glassfibre isolé avec Glassfibre dans l'ensemble, vernis en silicone partout. Bonne résistance a la température mais ne convient pas lorsque les fluides sont présents.

Spécifications**Specifications**

Product Code	WJ-033-D
General Description	Good temperature range but will not prevent ingress of fluids. Fairly flexible but does not provide good mechanical protection.
Cable Length	Priced per metre
Cable Type	Flat pair cable
Core / Strands	1/0.508mm
Outer Colour	Black
+ Colour	Yellow
- Colour	Blue
Insulation Material	Glassfibre insulated
Max. Temperature	+350°C (short periods up to 400°C)
Min. Temperature	-60°C
+ Conductor	Iron
- Conductor	Constantan



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Standards Met	BS
Thermocouple Type	J
Tolerance	Class 2