


Cable Thermocouple / Wire Type K PVC Insulated Flat Pair (BS)

Type K, 13/0.2mm
PVC Flat Pair



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Câble Thermocouple / Wire Type K PVC Insulated Flat Paire (BS)

L'image est pour l'illustration de codage de couleur seulement

Le prix est par metre - s'il vous plaît ajouter la quantité nécessaire pour le panier. Il ne reste plus que 613 metres en stock

Le câble de compensation ressemble a la caractéristique thermoélectrique du thermocouple pertinent sur une plage de température ambiante limitée, de 0 a 80 oC en général. Étant donné que ces câbles sont fabriqués a partir de matériaux a faible cout, des économies peuvent etre réalisées sur les installations de l'usine par rapport a l'exécution du véritable câble d'extension de thermocouple.

Spécifications**Specifications**

Product Code	WK-010-D
General Description	Good general purpose insulation for 'light' environments. Waterproof and very flexible.
Cable Length	Priced per metre
Cable Type	Flat pair construction
Core / Strands	13/0.2mm
Outer Colour	Red
+ Colour	Brown
- Colour	Blue
Insulation Material	PVC Insulated
Max. Temperature	+105°C
Min. Temperature	-10°C
+ Conductor	Nickel Chromium



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

- Conductor	Nickel Aluminium
Standards Met	BS
Thermocouple Type	K
Tolerance	Class 2