




LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

3 noyaux PFA isols, cran en cuivre tam, cble / fil de capteur PRT



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

3 noyaux PFA isols, cran en cuivre tam, cble / fil de capteur PRT

Haute temperature jusqu 260C, construction en tflon PFA extrud. Couramment utilis pour les thermomtres Pt100, convient galement dautres types tels que Pt130, Pt500 et Pt1000

Spcifications**Specifications**

Product Code	XE-2735
General Description	PFA insulated 7/0.2mm tin plated copper conductors, twisted cores, tin plated copper braided screen with PFA insulated outer jacket
Cable Length	10, 25, 50, 100m
Cable Type	PFA Teflon® construction
Core / Strands	7/0.2mm
Insulation Material	PFA Insulated
Max. Temperature	+260°C
Min. Temperature	-75°C
Number of Wires	3
Standards Met	Cores colour coded Red & White in accordance with IEC-751, black outer jacket