

Connecteur thermocouple miniature Prise en ligne FMTC-CU-FS Type Cu



Miniature in-line Socket
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Connecteur de thermocouple miniature Prise en ligne FMTC-CU-FS Type CU

Ideal pour une utilisation gnrale, tous les contacts sont polariss pour assurer une connexion correcte. Ces connecteurs acceptent des cbles thermocouples jusqu 4mm de diamtre.

Spifications**Specifications**

Product Code	XE-1372-001
General Description	Ideal for general purpose use. Contacts are polarised to ensure correct connection.
Thermocouple Type	IEC, ANSI, JIS Type CU
Connector Type	Socket
Connector Size	Miniature
Colour	White
Max. Temperature	220°C
+Positive Contact	Copper
-Negative Contact	Copper
Standards Met	BS EN 50212:1997. RoHS