


Connecteur thermocouple miniature Prise en ligne FMTC-CU-FS Type Cu



Miniature in-line Socket
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Connecteur de thermocouple miniature Prise en ligne FMTC-CU-FS Type CU

Ideal pour une utilisation gnrale, tous les contacts sont polariss pour assurer une connexion correcte. Ces connecteurs acceptent des cbles thermocouples jusqu 4mm de diamtre.

Spcifications**Specifications**

| | |
|----------------------------|---|
| Product Code | XE-1372-001 |
| General Description | Ideal for general purpose use. Contacts are polarised to ensure correct connection. |
| Thermocouple Type | IEC, ANSI, JIS Type CU |
| Connector Type | Socket |
| Connector Size | Miniature |
| Colour | White |
| Max. Temperature | 220°C |
| +Positive Contact | Copper |
| -Negative Contact | Copper |
| Standards Met | BS EN 50212:1997. RoHS |