


# Bayonet rÃfÃ©glable Pt100, Glassfibre Steel Inoxydable CAfÃ©ble surbradÃfÃ© - Type RTD / PRT



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Ajustable Bayonet PRT.**

Bayonet réglable (Pt100) PRT / RTD avec ressort 200mm avec bouchon de fixation réglable (11.5mm I.D. fente simple)

Disponible en longueurs 1M ou 2M.

**Specifications**

<b>General Description</b>	PRT with adjustable insertion depth
<b>Sensor Type</b>	Pt100 4 wire
<b>Cable Length</b>	1 or 2 metres
<b>Cable Type</b>	Glassfibre insulated lead with stainless steel Overbraid
<b>Core / Strands</b>	4 core 7/0.2mm
<b>Termination Type</b>	Bare Tails
<b>Max. Temperature</b>	+350°C
<b>Min. Temperature</b>	-50°C
<b>Standards Met</b>	IEC-60751 Class B
<b>Bayonet cap</b>	11.5mm I.D. single slot
<b>Spring</b>	6.0mm diameter x 200mm long