


Bayonet rÃfÃ©glable Pt100, Glassfibre Steel Inoxydable CAfÃ©ble surbradÃfÃ© - Type RTD / PRT



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Ajustable Bayonet PRT.

Bayonet réglable (Pt100) PRT / RTD avec ressort 200mm avec bouchon de fixation réglable (11.5mm I.D. fente simple)

Disponible en longueurs 1M ou 2M.

Specifications

General Description	PRT with adjustable insertion depth
Sensor Type	Pt100 4 wire
Cable Length	1 or 2 metres
Cable Type	Glassfibre insulated lead with stainless steel Overbraid
Core / Strands	4 core 7/0.2mm
Termination Type	Bare Tails
Max. Temperature	+350°C
Min. Temperature	-50°C
Standards Met	IEC-60751 Class B
Bayonet cap	11.5mm I.D. single slot
Spring	6.0mm diameter x 200mm long