


# Miniature Quick Wire Connecteur Thermocouple Socket FMTC-CU-FQ Type Cu



**Miniature in-line Socket - Quick Wire**  
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Miniature Quick Wire Connecteur Thermocouple Socket FMTC-CU-FQ Type CU**

Connexion jab-in facile, il suffit de foncer le fil et de serrer la vis. Les contacts sont polarisés pour éviter une polarisation incorrecte.

**Specifications****Specifications**

<b>Product Code</b>	XE-1382-001
<b>General Description</b>	Easy jab-in connection, just push in wire & tighten screw.
<b>Thermocouple Type</b>	IEC, ANSI, JIS Type CU
<b>Connector Type</b>	Quick Wire Socket
<b>Connector Size</b>	Miniature
<b>Colour</b>	White
<b>Max. Temperature</b>	220°C
<b>+Positive Contact</b>	Copper
<b>-Negative Contact</b>	Copper
<b>Standards Met</b>	BS EN 50212:1997. RoHS