


# Miniature Rectangulaire Thermocouple Connecteur Fascia Socket IM-N-FF Type N IEC



**Miniature Fascia Socket (with nylon mounting clip)**  
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Miniature Rectangulaire Thermocouple Connecteur Fascia Socket IM-N-FF Type N IEC**

Clip en nylon pour une fixation simple et rapide sur le panneau. Toutes les prises de montage de carnage sont polarisées au moyen de broches de taille inégale

**Spcifications****Specifications**

<b>Product Code</b>	XE-1257-001
<b>General Description</b>	All fascia mounting sockets are polarised by means of unequally sized pins
<b>Thermocouple Type</b>	N IEC
<b>Connector Type</b>	Fascia mounting socket
<b>Connector Size</b>	Miniature
<b>Colour</b>	Pink
<b>Max. Temperature</b>	120°C
<b>+Positive Contact</b>	Nickel / Chromium / Silicon
<b>-Negative Contact</b>	Nickel Silicon
<b>Standards Met</b>	BS EN 50212:1997. RoHS