

Thermomtre rsistance Pt100 Duplex 2 x 3 fils de classe A, tte KNE - Longueur de la sonde 100 mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications

Specifications

Product Code	XE-3655-001
Sensor Type	2 x Pt100 (100 Ohms @ 0°C) to IEC 751, Class A
Probe Diameter	6mm
Probe Length	100mm
Sheath Material	316 stainless steel
Head Type	KNE
Termination Type	6 way connection block for KNE head
Max. Temperature	+250°C (connection head @ 170°C)
Min. Temperature	-50°C
IP Rating	IP67
Accuracy	Class A
Number of Wires	2 x 3 Wire