



# LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

## Thermocouple isolante minérale de style industriel avec tête terminale KNE - Type K Inconel (6mm x 400mm)



6.0 X 100 mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

## Specifications

<b>Product Code</b>	XE-3248-001
<b>General Description</b>	Highly flexible, sheath can be bent/formed to suit many applications and processes
<b>Sensor Type</b>	Mineral insulated Thermocouple fitted with KNE terminal head
<b>Thermocouple Type</b>	Type K
<b>Thermocouple Junction</b>	Insulated hot junction
<b>Probe Diameter</b>	6.0mm diameter
<b>Probe Length</b>	400mm Length
<b>Sheath Material</b>	Inconel
<b>Head Type</b>	KNE, Cast aluminium finish, screw type lid with chain
<b>Gland</b>	M16 x 1.5mm Plated brass gland fitted
<b>Max. Temperature</b>	+1100°C
<b>Min. Temperature</b>	-40°C
<b>IP Rating</b>	IP67