

## Type K 1/0.711 Cble / fil thermocouple en fibre de verre double torsad en fibre de verre (ANSI)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Le prix est pour une longueur de 100 m

### Specifications

## Specifications

<b>Product Code</b>	XE-2812-100
<b>General Description</b>	Cores Twisted Together, 1/0.711mm Conductor Wire Size
<b>Cable Length</b>	100 metres
<b>Cable Type</b>	Fibreglass insulated twisted pair - High temperature
<b>Core / Strands</b>	2
<b>Conductor Cross Section (mm<sup>2</sup>)</b>	0.397
<b>Outer Colour</b>	n/a
<b>+ Colour</b>	Yellow
<b>- Colour</b>	Red
<b>Insulation Material</b>	Fibreglass - High Temperature
<b>Loop Resistance (Ohms per combined metre)</b>	2.5
<b>Max. Temperature</b>	+700°C
<b>Min. Temperature</b>	0°C
<b>Number of Wires</b>	3
<b>Overall Jacket Insulation</b>	No
<b>Overbraiding</b>	No
<b>Screen</b>	No



# LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

<b>Standards Met</b>	ANSI MC96.1 / wire tolerance: IEC 60584:1
<b>Thermocouple Type</b>	K (Nickel Chromium/Nickel Aluminium)
<b>Tolerance</b>	Class 1