


## Connecteur thermocouple rectangulaire miniature Prise de fascia FMTC-CU-FF Type Cu



**Miniature Fascia Socket (with nylon mounting clip)**  
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Connecteur de thermocouple rectangulaire miniature Prise de fascia FMTC-CU-FF Type CU**

Clip en nylon pour une fixation simple et rapide sur le panneau. Toutes les prises de montage de carnage sont polarises au moyen de broches de taille ingale

**Spcifications****Specifications**

<b>Product Code</b>	XE-1377-001
<b>General Description</b>	All fascia mounting sockets are polarised by means of unequally sized pins
<b>Thermocouple Type</b>	IEC, ANSI, JIS Type CU
<b>Connector Type</b>	Fascia mounting socket
<b>Connector Size</b>	Miniature
<b>Colour</b>	White
<b>Max. Temperature</b>	120°C
<b>+Positive Contact</b>	Copper
<b>-Negative Contact</b>	Copper
<b>Standards Met</b>	BS EN 50212:1997. RoHS