


# Plug thermocouple miniature quick Wire Connector Thermocouple Plug - Socket JM-T-MQ-FQ Type T JIS



**Miniature in-line Plug - Quick Wire**  
(dimensions in mm)



**Miniature in-line Socket - Quick Wire**  
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Miniature Thermocouple Connecteur Quick Wire Plug - Socket JM-T-MQ-FQ Type T JIS**

Connecteur miniature avec broches a bouchon plat. Connexion jab-in facile pour la résiliation rapide. Il suffit de pousser le fil et serrer la vis.

**Spécifications****Specifications**

<b>Product Code</b>	XE-1037-001
<b>General Description</b>	Easy jab-in connection for rapid termination
<b>Thermocouple Type</b>	T JIS
<b>Connector Size</b>	Miniature
<b>Connector Type</b>	Quick Wire Plug & Socket
<b>Colour</b>	Brown
<b>Max. Temperature</b>	220°C
<b>Pins</b>	Flat Pins
<b>+Positive Contact</b>	Copper
<b>-Negative Contact</b>	Copper Nickel
<b>Standards Met</b>	BS EN 50212:1997. RoHS