



# LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

## Thermocouple isolation minrale rponse rapide de 0,25 mm de diamtre avec prise miniature ANSI - Type K



IMAGE  
COMING  
SOON

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Veillez utiliser nos [formulaire de capteurs personnalis](#) pour soumettre vos besoins pour différentes longueurs dIM, car tous terminés par une prise une connexion simple et rapide.

### **Vous avez besoin de calibrer votre capteur de température ?**

Traable selon les normes dtalonnage ISO17025 (UKAS), Labfacility offre un choix dtalonnage traable en 3 ou 5 points laide de notre installation dtalonnage interne.

**Veillez noter que si vous ajoutez une option de calibrage votre commande, vous serez contact pour les informations de livraison.**

### **Spifications**

## **Specifications**

<b>General Description</b>	Highly flexible, sheath can be bent/formed to suit many applications and processes
<b>Sensor Type</b>	Type K ANSI MB STYLE probes 0.25mm
<b>Thermocouple Junction</b>	Single element, junction insulated or grounded
<b>Probe Length</b>	100mm / 150mm / 250mm
<b>Sheath Material</b>	Stainless steel 316 sheath
<b>Termination Type</b>	Fitted Miniature Plug Type ANSI
<b>Range</b>	-40 to +260°C
<b>Important Information</b>	A 1%+ Length Tolerance applies