


# PCB Montage Thermocouple Connector Socket IM-K-PCB Type K IEC



**Miniature PCB Socket**  
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**PCB Montage Thermocouple Connecteur Socket IM-K-PCB Type K IEC**

Spécialement conçu pour les circuits imprimés, les broches solides pour le montage direct sur carte et la disposition pour capteur CJC

**Specifications****Specifications**

<b>Product Code</b>	XE-1109-001
<b>General Description</b>	Versatile, high quality circuit board mounting miniature socket
<b>Thermocouple Type</b>	K IEC
<b>Connector Type</b>	PCB Socket
<b>Connector Size</b>	Miniature
<b>Colour</b>	Green
<b>Max. Temperature</b>	220°C
<b>Pins</b>	Brass
<b>+Positive Contact</b>	Nickel Chromium
<b>-Negative Contact</b>	Nickel Alloy
<b>Standards Met</b>	BS EN 50212:1997. RoHS