

## Lascar EL-P-TPX+ Calibratable thermistor probe



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

High accuracy calibratable, 3 metre long thermistor probe with stainless steel tip for use with EL-WiFi-TPX+ ([XE-4076-001](#))

### Specifications

## Specifications

<b>Product Code</b>	XE-4081-001
<b>Model</b>	EL-P-TPX+
<b>Type</b>	Calibratable thermistor probe
<b>Accuracy</b>	+/-0.2°C (-15 to +80°C) / +/-0.6°C (-40 to 125°C)
<b>Cable Type</b>	3m (118.1")
<b>Measuring range</b>	-40 to +125°C / -40 to + 257°F