


PCB Montage Thermocouple Connector Socket IM-T-PCB Type T IEC



Miniature PCB Socket
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

PCB Montage Thermocouple Connecteur Socket IM-T-PCB Type T IEC

Spécialement conçu pour les circuits imprimés, les broches solides pour le montage direct sur carte et la disposition pour capteur CJC

Specifications**Specifications**

Product Code	XE-1199-001
General Description	Versatile, high quality circuit board mounting miniature socket
Thermocouple Type	T IEC
Connector Type	PCB socket
Connector Size	Miniature
Colour	Brown
Max. Temperature	220°C
+Positive Contact	Copper
-Negative Contact	Copper Nickel
Standards Met	BS EN 50212:1997. RoHS