



# LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

**Diamtre de 0,75 mm avec queues nues - isol ou  
mis la terre**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Sondes thermocouples joint de transition de diamtre fin**

Les jonctions de thermocouple sont disponibles en version mise la terre, pour un temps de rponse plus rapide et un cot infrieur ; ou isol, pour lisolation lectrique du thermocouple de la gaine.

Si vous avez besoin de plus grandes quantits, veuillez sales@labfacility.com envoyer un e-mail ou tlphoner notre ligne de vente 01243 871280.

**Vous avez besoin de calibrer votre capteur de temprature ?**

Traable selon les normes dtalonnage ISO17025 (UKAS), Labfacility offre un choix dtalonnage traable en 3 ou 5 points laide de notre installation dtalonnage interne.

**Veillez noter que si vous ajoutez une option de calibrage votre commande, vous serez contact pour les informations de livraison. .**

**Spfications****Specifications**

<b>Sensor Type</b>	Mineral Insulated
<b>Thermocouple Type</b>	K
<b>Element/hot junction</b>	Single element, junction insulated or grounded from sheath
<b>Probe Diameter</b>	0.75mm
<b>Probe Length</b>	100,150, 250mm
<b>Sheath Material</b>	K - Stainless steel 316
<b>Termination Type</b>	Fortfix sealed with = +VE leg 15mm long and -VE leg 10mm long
<b>Max. Temperature</b>	+650°C
<b>Min. Temperature</b>	-40°C
<b>Tolerance</b>	Class 1
<b>Standards Met</b>	In accordance with IEC 584