

Minuterie de thermomtre de cuisson numrique avec alarme et aimant arrire



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

TypeT-1/0.2

Specifications

Product Code	XE-4350-001
Model	LDT-100
General Description	Digital Cooking Thermometer Timer With Alarm & rear magnet
Type	Timer-Thermometer 2 in 1
Accuracy	+/-1°C (0-10°C) /-1.5°C(101-25°C)
Readings	Timer range: 99 Min 59 Sec
Resolution	0.1°C(0.2°F)
High Alarm Level	High temperature alarm function
Battery	1 x 1.5V(AAA)
Cable Length	Wire: 120cm length
Length	10m
Colour	Metal finish
Dimensions	6.5 x 7 x 1.7cm
Functions	°C / °F switch
Temperature Range	0-250°C / 32-482°F
Probe Material	304 stainless steel
Weight	100g
Construction	Eco-friendly ABS plastic + Stainless steel front shell
Important Information	Hanging & flip stand