

## Minuterie de thermomtre de cuisson numrique avec alarme et aimant arrire



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

TypeT-1/0.2

## Specifications

<b>Product Code</b>	XE-4350-001
<b>Model</b>	LDT-100
<b>General Description</b>	Digital Cooking Thermometer Timer With Alarm & rear magnet
<b>Type</b>	Timer-Thermometer 2 in 1
<b>Accuracy</b>	+/-1°C (0-10°C) /-1.5°C(101-25°C)
<b>Readings</b>	Timer range: 99 Min 59 Sec
<b>Resolution</b>	0.1°C(0.2°F)
<b>High Alarm Level</b>	High temperature alarm function
<b>Battery</b>	1 x 1.5V(AAA)
<b>Cable Length</b>	Wire: 120cm length
<b>Length</b>	10m
<b>Colour</b>	Metal finish
<b>Dimensions</b>	6.5 x 7 x 1.7cm
<b>Functions</b>	°C / °F switch
<b>Temperature Range</b>	0-250°C / 32-482°F
<b>Probe Material</b>	304 stainless steel
<b>Weight</b>	100g
<b>Construction</b>	Eco-friendly ABS plastic + Stainless steel front shell
<b>Important Information</b>	Hanging & flip stand