



# LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

## JIS Type K Thermocouple pointe soude expose 1 / 0.3mm PFA Twin Twisted Cable avec queues nues 10m



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Fin de la queue nue.

Jonction soude et expose pour des temps de rponse rapides. Construction durable avec isolation PFA pour une rsistance aux huiles, aux acides, aux fluides et une bonne rsistance mcanique et flexibilit.

## Specifications

|                            |   |
|----------------------------|---|
| <b>Product Code</b>        | XE-0527-001   |
| <b>General Description</b> | Welded, exposed junction for fast response times. Durable construction with PFA insulation for resistance to oils, acid, fluids and good mechanical strength and flexibility. |
| <b>Manufacturing Code</b>  | Z2-K-10.0-1/0.3-T-JIS   |
| <b>Sensor Type</b>         | Exposed Junction  |
| <b>Thermocouple Type</b>   | JIS Type K  |
| <b>Cable Length</b>        | 10m (32.8ft)  |
| <b>Cable Type</b>          | PFA Twin Twisted  |
| <b>Core / Strands</b>      | 1/0.3mm   |
| <b>Cable Termination</b>   | Bare Tails  |
| <b>Max. Temperature</b>    | +250°C  |
| <b>Min. Temperature</b>    | -75°C   |
| <b>Series</b>              | Exposed Junctions   |