



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

JIS Type K Thermocouple pointe soude expose 7 / 0.2mm PFA Cble torsad double avec queues nues 2m



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Fin de la queue nue.

Jonction soude et expose pour des temps de rponse rapides. Construction durable avec isolation PFA pour une rsistance aux huiles, aux acides, aux fluides et une bonne rsistance mcanique et flexibilit.

Specifications

Product Code	XE-0534-001
General Description	Welded, exposed junction for fast response times. Durable construction with PFA insulation for resistance to oils, acid, fluids and good mechanical strength and flexibility.
Manufacturing Code	Z2-K-2.0-7/0.2-T-JIS
Sensor Type	Exposed Junction
Thermocouple Type	JIS Type K
Cable Length	2m (6.56ft)
Cable Type	PFA Twin Twisted
Core / Strands	7/0.2mm
Cable Termination	Bare Tails
Max. Temperature	+250°C
Min. Temperature	-75°C
Series	Exposed Junctions