




# LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

## IEC Type K 11.8kg Pull Power (fer cheval) Thermocouple magnétique, cble isol PFA avec tresse en acier inoxydable se terminant en miniature ou fiche standard



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

## Specifications

<b>General Description</b>	11.8kg Pull Power (Horseshoe) Magnet Thermocouple
<b>Sensor Type</b>	Magnet Thermocouple
<b>Thermocouple Type</b>	Type K
<b>Magnet</b>	45 x 30 x 30mm (L x H x W) - Max pull: 11.8 KG
<b>Cable Length</b>	1.5M or 3.0M
<b>Cable Type</b>	PFA Insulated Cable with Stainless Steel Over-Braid
<b>Core / Strands</b>	7/0.2mm
<b>Cable Termination</b>	Miniature or Standard Plug
<b>Max. Temperature</b>	250°C
<b>Min. Temperature</b>	-50°C
<b>Standards Met</b>	IEC