




LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

IEC Type K 11.8kg Pull Power (fer cheval) Thermocouple magnétique, cble isol PFA avec tresse en acier inoxydable se terminant en miniature ou fiche standard



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications

General Description	11.8kg Pull Power (Horseshoe) Magnet Thermocouple
Sensor Type	Magnet Thermocouple
Thermocouple Type	Type K
Magnet	45 x 30 x 30mm (L x H x W) - Max pull: 11.8 KG
Cable Length	1.5M or 3.0M
Cable Type	PFA Insulated Cable with Stainless Steel Over-Braid
Core / Strands	7/0.2mm
Cable Termination	Miniature or Standard Plug
Max. Temperature	250°C
Min. Temperature	-50°C
Standards Met	IEC