



**Thermocouple magnétique bouton de traction IEC  
Type K de 1,9 kg, cble isol PFA avec tresse en  
acier inoxydable se terminant en miniature ou en  
prise standard**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

## Specifications

<b>General Description</b>	1.9kg Pull Button Magnet Thermocouple
<b>Sensor Type</b>	Magnet Thermocouple
<b>Thermocouple Type</b>	Type K
<b>Magnet</b>	19.1mm Diameter / 12.7mm height. Max Pull: 1.9kg
<b>Cable Length</b>	1.0m or 2.0m
<b>Cable Type</b>	PFA Flat Pair with stainless steel overbraid
<b>Core / Strands</b>	7/0.2mm
<b>Termination Type</b>	Miniature or Standard Plug
<b>Max. Temperature</b>	250°C
<b>Min. Temperature</b>	-50°C
<b>Standards Met</b>	IEC