



**Thermocouple magnétique IEC Type K 4,5 kg (fer cheval), cble isol PFA avec sur-tresse en acier inoxydable se terminant par une prise miniature**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

## Specifications

<b>General Description</b>	4.5kg Pull Power (Horseshoe) Magnet Thermocouple
<b>Sensor Type</b>	Magnet Thermocouple
<b>Thermocouple Type</b>	Type K
<b>Magnet</b>	30 x 20 x 20mm (L x H x W). Max Pull: 4.5kg
<b>Cable Length</b>	1.0M or 2.0M
<b>Cable Type</b>	PFA insulated lead with stainless steel overbraid
<b>Core / Strands</b>	7/0.2mm
<b>Termination Type</b>	Miniature or Standard Plug
<b>Max. Temperature</b>	250°C
<b>Min. Temperature</b>	-50°C
<b>Standards Met</b>	IEC