



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Patch en caoutchouc de silicone thermocouple PFA twin torsad avec terminaison queues nues - Type K JIS - 2m



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications**Specifications**

Product Code	XE-0559-001
General Description	Resistant to many oils & chemicals
Sensor Type	Thermocouple junction encapsulated in a silicone rubber patch
Thermocouple Type	K JIS
Cable Length	2m
Cable Type	PFA insulated twin-twist wire 1/0.2mm
Termination Type	Bare Tails
Max. Temperature	+150°C
Min. Temperature	-50°C
Standards Met	JIS (Japanese Industrial Standards)
Colour	Positive: Red / Negative: White
Dimensions	Length 40mm x Width 13mm x Height 5mm Rubber patch