



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

ANSI Type K 0.7kg Pull Button Magnet Thermocouple, PFA isol avec prise miniature



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications**Specifications**

General Description	Ideal for ferrous metal surface measurement
Thermocouple Type	K
Thermocouple Junction	Grounded at tip
Magnet	1/2" diameter button type, 0.7kg Pull
Cable Length	0.5m, 1m, 2m
Cable Type	PFA insulated twin-twist wires
Cable Termination	Miniature Fitted Plug
Max. Temperature	+250°C (plug 200°C)
Standards Met	ANSI MC96.1
Colour	+ Red / - Yellow