

Capteurs de temprature de plaquette de frein thermocouple JIS - Type K



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Ces capteurs de thermocouple doivent être utilisés pour surveiller la température des plaquettes de frein jusqu'à +400°C.
Ils sont conçus pour être intégrés dans les pads.

Chaque capteur a 1 m de câble isolé en fibre de verre, avec Overbraid en acier inoxydable. Disponible en thermocouple de type K et avec un choix de diamètres de pointe. Chaque pointe est scellée avec une pâte céramique haute température (Fortafix) qui peut résister à des températures allant jusqu'à +400°C.

Specifications

Specifications

| | |
|------------------------------|---|
| General Description | Designed for monitoring the temperature of brake pads |
| Sensor Type | Grounded Junction |
| Thermocouple Type | K |
| Probe Diameter | 3.2mm, 4.0mm or 6.0mm |
| Cable Length | 1m |
| Cable Type | Fibre glass with stainless steel Overbraid 7/0.2mm |
| Termination Type | Bare Tails |
| Max. Temperature | +400°C |
| Standards Met | JIS (Japanese Industrial Standards) |
| Colour | Positive: Red / Negative: White |
| Important Information | Potting sealant: Ceramic paste (Fortafix) |